MATRIX SURGICAL USA

Where Art Meets Science

CE 0120

STERILE

Sterility guaranteed unless package is opened, damaged or wet.

Distributed by
Matrix Surgical USA
4025 Welcome All Road, Suite 120,
Atlanta, GA 30349 USA
Tel: 1-404-855-4592 Fax: 1-404-344-0283
www.matrixsurgicalusa.com

Jewel Management Ltd.
Crown House, Waller Road, Devizes, Wiltshire
SN10 2GH, United Kingdom

CONFORMER IFU

FOR SINGLE USE ONLY

PRODUCT INFORMATION - English

DESCRIPTION

Matrix Surgical USA Ocular Conformers are manufactured of acrylic and maintains the shape of the eye and prevent closure or adhesion during the healing process. The Ocular Conformers are sold vented or non-vented.

INDICATIONS

Conformers are temporary devices and intended as post-operative enucleation stints to provide a smooth surface for the lids to blink over while protecting the surgical closure beneath. When healing is complete, conformers are replaced with prosthesis. Conformers can also be used for eviscerations and symblepharon repair.

CONTRAINDICATIONS

The following are specific contraindications, adverse effects, and warnings that must be understood by the surgeon and explained to the patient. The general surgical risks that are not included should also be explained to the patient prior to surgery.

1. Active infections.
2. Patients who are unwilling or unable to follow postoperative instructions due to accompanying conditions (mental or physical).
3. Foreign body sensitivity. Where material sensitivity is suspected, appropriate tests should be made to rule out this possibility prior to implantation.
4. Use in tissue that has been compromised by cancer therapies.
5. Limitations in blood supply and/or systemic disorders that may slow healing and increase the possibility of infection and/or rejection of the conformers.
6. Any degenerative disease process that would adversely affect proper placement of the conformers.
7. Inadequate coverage with healthy tissue.

POSSIBLE ADVERSE EFFECTS

1. Vascular changes.
2. Allergic reaction to the conformer.
3. Dehiscence of the conjunctiva or tenons.
4. Extrusion of the conformer.
5. Tissue adhesion to the conformer.
6. Pain, discomfort, and/or abnormal sensation due to presence of the conformer.
7. Superficial and/or deep infection.

WARNINGS

1. Conformers are never to be re-used or re-sterilized.
2. Adequate patient instruction is a key factor in determining the success of the surgical procedure. Postoperative follow-up and care are very important.
3. The surgeon must be thoroughly familiar with the conformers, the method of application, the instrumentation, and the surgical procedure. The patient must be made aware of the risks of using the conformers including possible adverse effects.
4. Responsibility for completion of adequate training, proper patient selection, and choice and placement of the conformers rests with the physician.

INSTRUCTIONS FOR USE

The following is a guide for informational uses only. Since a number of surgical techniques may be employed, the surgeon may wish to alter details in accordance with his or her clinical experience and medical judgment. The selection of the conformer and its composition, size, and shape is made by the surgeon after evaluating the anatomical and functional requirements of the patient.

During the closure of an evisceration or enucleation procedure, place the conformer into the socket after the conjunctiva layer has been sutured and closed.

Remove the conformer after edema of the socket has subsided and proceed with Post-Operative procedures.

WARRANTY

All products are warranted to be free from defects in material and workmanship. No warranty is made for any purpose other than in the product specifications and labeling.

CAUTION

Federal law (USA) restricts this device to sale by or on the order of a physician.

Additional information pertaining to specific devices or issues may be obtained from Matrix Surgical USA.

© 2015 Matrix Surgical USA.

The MATRIX SURGICAL service mark is owned by or used under the license and authority of Matrix Surgical Holdings, LLC.

Labeling Symbols

	Attention, See Instructions for Use
	Do Not Reuse
	Sterilized using Ethylene Oxide
	Use by Date
	Legal Manufacturer
	Authorized Representative in the European Community
	Batch Code / Lot Number
	Keep Dry
	Do Not Use if Package is Damaged

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